

TECHNICAL DATA SHEET

DETER BIANCO

CHARACTERISTICS

Optical brightener for all types of tissues

The product consists of white granules

APPLICATIONS AND DOSAGE

Dissolve from 3 to 5 grams approximately of DETER BIANCO per liter of warm water.

The product is already active at 30°C and it increases its effectiveness with the gradual increasing of temperature.

The ideal temperature is from 40/50 ° C.

Once mixed with water, leave the items into it.

Leave for about 20 minutes, move often the items for a uniform action.

At the end of treatment rinse thoroughly.

CONTENTS/TECHNICAL DATA

Appearance and colour: white powder

Odour: odourless and slightly pungent

pH (1% in H₂O): 7,5 to - 10

DIRECTIONS

Store in a dry place and keep container tightly closed.

Le nostre circolari servono esclusivamente a dare informazioni generali sull'uso dei nostri prodotti. Si consiglia di vagliare responsabilmente qualsiasi influenza possa verificarsi sui capi. Le ns. indicazioni riflettono le ns. cognizioni tecniche più aggiornate, ciò non implica alcuna garanzia legale da parte nostra. Caratteristiche e dati possono subire variazioni senza preavviso da parte nostra.

G.B.M. Elettrochimica S.R.L.

Sede Legale – Via Fiumicino San Mauro 120/130 – 47039 Savignano sul Rubicone (FC)

Tel. 0541 930058 r.a. – Fax Comm. 0541 810218 - Fax Amm.ne 0541 815280 - www.gbmprodottichimici.it - gbm@gbmprodottichimici.it

Capitale Sociale € 20.000,00 i.v. - Cod. Fisc. e P.Iva 00958710402 – R.E.A. 192182